

PRODUCT PROFILE

NAME OF PRODUCT: Vitallium® 2000 Partial Denture Frameworks

MANUFACTURER: Dentsply Austenal

BRIEF DESCRIPTION: Premium Chromium Cobalt Alloy


BENEFITS FOR DOCTOR'S OFFICE: Vitallium's 2000's tensile strength of over 855 MPa's makes it extremely fracture-resistant and with an elongation of 9%, the partials adjust like gold.

BENEFITS FOR PATIENT: Strength, superior fracture resistance, and smaller, lighter design

WHAT TO SEND LABORATORY: Full arch impression after tooth modification, opposing and registration. Rx should include framework design.

INDICATIONS: All removable partial denture frameworks

CONTRA-INDICATIONS: None (are Beryllium-free)



- PREPARATION INSTRUCTIONS: May send in study model for survey and design recommendations.
Prepare rest seats and guide plans
- DAYS IN LAB: 4
- SEATING INSTRUCTIONS: Adjustments should be none to minimal.
- AFTER CARE: Remove and brush appliance after meals and at bedtime.
- INSURANCE CODE: D5213 Maxillary Partial Denture - cast metal frame with resin denture bases (including any clasps, rests, teeth)
D5214 Mandibular Partial Denture - cast metal frame with resin denture bases (including any clasps, rests, teeth)
- CHAIRSIDE INSTRUCTIONS: Routine